


Biocidal ZF™

Spray Disinfectant for Incubators and Sterile Cabinets in Cell Culture Area

Cell Culture Area



Effectiveness:

Biocidal ZF™ is active against bacteria, fungi, spores and enveloped viruses.

Non-volatile:

Active ingredients are non-volatile. They protect cell cultures from microbial contamination and do not invade cell cultures via air. Thus cell cultures are protected against contamination and the disinfectant itself.

Application for incubators

- Empty incubator
- Spray Biocidal ZF™ on all surfaces and let settle
- Put cultures right back

Surfaces will be protected for 2-4 weeks. (Before disinfecting you may clean incubator with Biocidal ZF™ if necessary)

Application in sterile benches

- Empty sterile bench
- Then clean with Biocidal ZF™ and finally spray on all surfaces.

Protection will prevail until next workshift. No combustible gases can form, unlike by use of ethanol.

Classification according to EU-Regulation
No. 1272/2008

Not hazardous

Active ingredients:

Quaternary ammonium compounds
-free of aldehydes-

Storage:

Room temperature

Materials compatibility:

Stainless steel

Cat.#

WAK-ZF-1	1 x 1	liter spray bottle
WAK-ZF-6	6 x 1	liter spray bottle
WAK-ZF-12	12 x 0.2	liter spray bottle

L-01-136-02-0-01

Siemensstraße 9
D-61449 Steinbach/Ts. • Germany
Tel. +49 61 71 / 28 43 0 • Fax +49 61 71 / 69 56 46
e-mail: info@wak-chemie.com
www.wak-chemie.com



**WAK - Chemie
Medical GmbH**

