


# Bandit™

STEM Electrophoresis Kit  
QP-1400-01

## TECHNICAL SPECIFICATIONS



<b>GEL SIZE</b>	55 x 60 mm
<b>GEL VOLUME</b>	30 ml
<b>BUFFER VOLUME</b>	30 ml
<b>COMB</b>	6 or 9 wells
<b>ELECTRODE</b>	Corrosion-resistant alloy
<b>INPUT VOLTAGE</b>	5V DC, USB-C
<b>OUTPUT VOLTAGE</b>	40V DC
<b>SAFETY FEATURES</b>	Low voltage, over-voltage detection, active current limiting
<b>WEIGHT</b>	6.4 oz

1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com